

Portable Autoclave 'P' Type/Cooker Type

- All body made of Aluminium
- Ring and Lid made of aluminium Casting,
- used for sterilization Surgical instruments and dressing.
- Material – Aluminium Hydraulically tested upto 30 psi.
- Fitted with pressure gauge
- The unit has spring loaded safety valve of pressure adjustment and a dead weight safety valve

INCLUDING:

Material	Aluminium
Orientation	Vertical
Pack Type	Box
Usage	Laboratory
Capacity	20 L, 25L



Technical Specification of Portable Autoclave 'P' Type/Cooker Type



Model	Size(mm)	Size(inches)	Capacity(liter)
901(A)	225x275	09x11	13ltr
901(A)	300x300	12x12	20ltr
901(A)	300x500	12x20	40ltr

CONTACT:



Plot No. 6 - Ekta Vihar, Ambala Cantt.133001



Mobile No.: 8295732110, 7056732110



Sales@labsolutionindia.com



Labsolutionindia.com



GeM
Government
e Marketplace
OEM Regd

