



RESEARCH MEDICAL **MICROSCOPE**

RESEARCH MEDICAL MICROSCOPE

- ✓ Frame– All aluminum pressure die casted body.
- ✓ Paint – P.U. coating in light cream shade.
- ✓ Nose piece– Quadruple nosepiece.
- ✓ Focusing system– Separate coarse ultra smooth fine motion.
- ✓ Mechanical stage– Mechanical stage for smooth functioning.
- ✓ Condenser – Condenser Bright field sub stage condenser having 1.2 N.A.
- ✓ Objectives– 10x, 40x, and 100x and a paired set of 10x wide field eyepieces.
- ✓ Eye piece– WF-10X

About Research Medical Microscope

standard configurations

Objectives 10x, 40x, and 100x
Eyepiece WF-10X

Optional accessories

- Wooden Box
- 2/3 Mega Pixel Digital Camera for Microscope
- Pin Type Lamp (Bulb/Led) 220V
- H10x/H15

