

AUTOCLAVE VERTICAL DELUXE



Description: The construction is as per vertical autoclave designed to meet the Quality of any international standard. The lid & locking system is designed in such a way to be handled by a square handle at top. All the locking system and lid is made of Thick Gun Metal duly chrome plated. The outer is made of square Shape. Water level & pressure is locked in an enclosure. The castor wheels are provided at the bottom. Supplied with complete IQ, OQ & PQ documentation.

CONSTRUCTION AND SALIENT FEATURES

- The outer body is rectangular in shape.
- The lid and tightening system is provided of gun metal duly chrome plated.
- It is complete with all the optional accessories thus having triple safety system.
- It is supplied with extra heater which can automatically be used as & when required with separate switch.
- Automatic water inlet valve.
- Temperature controller cum indicator.
- Validation port.
- We can match this model with any world class autoclave.

CAPACITY

- 22 Litres
- 40 Litres
- 50 Litres
- 78 Litres
- 98 Litres
- 150 Litres

OPTIONS

- Triple wall
- Double wall
- PLC controller instead of digital controller