

# Digital Type Oil Bath

Digital Type Oil Bath is double walled with inside make up of S.S. and outside of m.s. sheet duly enamel painted. The gap between the two walls is fitted with high grade glass wool to minimize thermal loss. Temperature range from  $50^{\circ}\text{C}$  to  $250^{\circ}\text{C} \pm 2^{\circ}\text{C}$  is controlled by "EGO"/JUMO" GERMAN hydraulic capillary type thermostat. Or digital temperature PID controller.

## About:

These oil baths can be availed with added accessories such as GMP model, and digital indicator and controller for constant temperature. Further, it is provided with a heavy duty motor that makes its operation swift and easy.



# Key Features Of Digital Type Oil Bath



Microprocessor Based PID Digital Temperature Indicator-cum-Controller maintains maximum temperature accuracy.

User's Self-compensation function to control difference between real temperature and indicated temperature.

Stainless Steel Bath for superior Durability & High Thermal efficiency.  
Flat Stainless Steel Lid keeps temperature constant.

Over temperature & Over current protection.

Industrial grade energy efficient heaters reduce power bills drastically.  
Double walled construction.

Space between the walls is fitted with special grade ceramic fibre in order to avoid thermal losses.

MOC (Inner) : Inner bath chamber and lid made of Stainless Steel (SS-304 grade).

MOC (Outer) : Outer chamber is made of Mild Steel duly powder coated  
To work on 220/230V AC supply.

## CONTACT:



Plot No. 6 - Ekta Vihar, Ambala Cantt.133001



Mobile No.: 8295732110, 7056732110



Sales@labsolutionindia.com



Labsolutionindia.com



**GeM**  
Government  
e Marketplace  
OEM Regd

