

Dry Bath INCUBATOR

The Digital Dry Baths offer digital control over both temperature and time, eliminating the need for an external thermometer or timer.

The cavities within the high-grade aluminum blocks are precision-machined to match the conical size and shape of various tube sizes, providing unsurpassed temperature uniformity.

ABOUT

- Block size:
5 ml x 24
- Power consumption :
65 watts
- Temperature variation:
+₋ 1c
- Temperature setting range:
37c +₋ 1c
{adjustable range 30c to 80c
- Temperature sensor:
Solid state electronic
- Temperature sensor :
For high accuracy and linearity
- Supply :
230v ac+ ₋ 5%, 50 hz



Key Features

Dry Bath Incubator



Body	Metal Body
Temp. Controller	Digital
Temperature Range	Ambient +5 to 100 degrees C
Accuracy	± .5 Degrees C
Heating Block	24 holes of 12.7 mm diameter
Timer	1 to 999 minutes
Heating Time	15 Minutes from RT to max.
W x D x H	168 (W) x 232 (D) x 93 (H)
Weight	2.5 Kgs.
Mains Input	230V/50Hz A.C. Supply

CONTACT:

-  Plot No. 6 - Ekta Vihar, Ambala Cantt.133001
-  Mobile No.: 8295732110, 7056732110
-  Sales@labsolutionindia.com
-  Labsolutionindia.com



GeM
Government
e Marketplace
OEM Regd

