

# Deluxe Blood Donor Couch

Deluxe Blood Couch has solid surface seat to ensure easy & hygienic premium rated foam, cushion and spring component for donor comfort. With the help of electric power, we can fully adjust the donor position in between head low-feet high to head high-feet lower position and arm rest (left and right) swing out for easy seating as well as up & down of chair such position minimization of donor.

## Specification:

- Dimension: Minimum 60 inch (L) x 36 inch (W)
- Head height reachable up to approx. 45 inch
- Foot height reachable up to approx. 35 inch
- Seat Height approx. 25 inch
- Capacity of patient weight. 150 kg
- Power: 220 V AC, 50 Hz, Single Phase.



# Key Features

## Of Deluxe Blood Donor Couch



Motorized Donor Chairs suitable to use for blood banks, infusion treatment centers, and hospital, clinics, and aphaeresis, plasma centers and other similar discipline.

Remote switch for easy repositioning.

All hardwood frame construction.

Push pillow back for donor comfort.

Lockable castors for easy mobility.

Fire Resistant Upholstery & Punishment proof laminates.

High quality contract grade vinyl upholstery with antimicrobial/antibacterial coating to prevent from stain, abrasion and bacterial growth.

All exposed metal ports are chrome plate for easy cleaning & durability.

## CONTACT:

 Plot No. 6 - Ekta Vihar, Ambala Cantt.133001

 Mobile No.: 8295732110, 7056732110

 Sales@labsolutionindia.com

 Labsolutionindia.com



**GeM**  
Government  
e Marketplace  
OEM Regd

