

# BENCH TOP SHAKER

These bench top shakers incubating orbital shakers are available with and without heating or refrigeration. The clear lid lifts up and out of the way providing visibility and giving easy access to all samples.

01 > 02 > 03 > 04



Large easy to read displays shows set point and actual values.

Integral timer, 0,1 to 99,9 h, shows elapsed time or counts down, automatic switch-off and audible alarm when program is complete.

Micro process or controller maintains precise temperature, stores program for repeat use and recalls last set points, even when unit has been shut off.

Air circulation system ensures rapid temperature equilibration and uniformity, multifunction reservoir humidifies the chamber to minimize sample evaporation.

# Key Features Of Bench Top Shaker



Orbit Dia/Tilt Angle	22 mm
Speed Range	30-350 RPM
Max. Load	2.5 KGs.
Time	1 minute to 20 hours
Speed display	LCD
Timer display	LCD
Power	20 W
Dimension of Platform	315x315 mm
Equipment Dimension	420x360x125 mm
Weight	7 Kgs.
Power Supply	220 volts AC

## CONTACT:

 Plot No. 6 - Ekta Vihar, Ambala Cantt.133001

 Mobile No.: 8295732110, 7056732110

 Sales@labsolutionindia.com

 Labsolutionindia.com



**GeM**  
Government  
e Marketplace  
Regd. OEM

