

# Blood Donor Couch

The Blood Donor couch is designed to provide a comfortable position to the donor during blood donation process. Blood Donor position can be shifted from upright body position in case of vasagonal attack. The Meditech Blood Donor Couch is a fully upholstered Blood donor chair.

## About:

Solid surface seat to ensure maximum comfort for the donor. The seat has a premium rated foam filled back rest supported by a steel structure.

The sturdy construction of this seating ensures that the chair can be used year after year.



# Key Features Of Blood Donor Couch



Electric blood donor Chair Fully Comfortable.

Body Material – Wooden , Steel , with Rexine Cover for long life use.

Fully Automatic Remote Function for backrest & foot rest .

220 VAC Operated 50 Hz.

Arm Rest two Side with Height Adjustable.

Long life Mechanical System with easy use & Protection Limits Switch.

Ergonomic designed.

Comfortable and right posture.

Light weight and stylish design.

## CONTACT:

 Plot No. 6 - Ekta Vihar, Ambala Cantt.133001

 Mobile No.: 8295732110, 7056732110

 Sales@labsolutionindia.com

 Labsolutionindia.com



**GeM**  
Government  
e Marketplace  
OEM Regd

