

# Blood Roller Mixer

This is an equipment to make blood homogeneous. It is shake-less and gentle mixing process. It is essential pre-process for perfect result of RBC, WBC and platelet counting on Hematology Analyzer.

## About:

- For thorough mixing of blood samples and cell suspensions.
- Keeps the contents of closed tubes in gentle motion by rocking while rolling.
- Compact and ergonomic design.
- 4 rolls of PVC rods for smooth operation.
- Can accommodate various sizes of centrifuge & tubes.



# Key Features

## Of Blood Roller Mixer



Rolling Speed	10 – 30 RPM
Rocking Frequency	5 Cycles per minute at 35° angle
No. of Rollers	4 or 6
Dimensions	455 x 265 x 120 mm
Exterior Body	MS Coated
Rolling Rods Material	PVC, each of 320 mm
Maximum Load	5 kg.
Supply	220 / 230 V, 50 Hz
Timer	No

### CONTACT:



Plot No. 6 - Ekta Vihar, Ambala Cantt.133001



Mobile No.: 8295732110, 7056732110



Sales@labsolutionindia.com



Labsolutionindia.com



**GeM**  
Government  
e Marketplace  
OEM Regd

