

Ice Lined Refrigerator

With the aim of delivering qualitative products, we are engaged in offering our clients with an excellent quality Ice Lined Refrigerator For Vaccine Storage. These are widely demanded in clinics, hospitals, blood banks, research institutes etc. and are ideal for storing vaccines.

About:

- We are offering our clients with an excellent quality Ice Lined Refrigerators For Vaccine Storage.
- These are widely demanded in clinics, hospitals, blood banks, research institutes etc. and are ideal for storing vaccines.



Key Features

Of Ice Lined Refrigerator



Temperature Range	2 to 8°C
Temperature Display	LCD
Temp. Control	Microprocessor Based
Temp. Holding Time	24 Hours from 2 to 8 degrees
Construction	Triple Walled
Outer Body	Mild Steel with epoxy powder coating
Inner Chamber	Stainless Steel
Baskets	SS wire mesh baskets
Insulation	75 mm thick puff insulation
Sub Door	Made of 8 mm thick acrylic
Compressor	Hermetically sealed with CFC free refrigerant
Supply	220 volts AC

CONTACT:

-  Plot No. 6 - Ekta Vihar, Ambala Cantt.133001
-  Mobile No.: 8295732110, 7056732110
-  Sales@labsolutionindia.com
-  Labsolutionindia.com



GeM
Government
e Marketplace
OEM Regd

