

# LABORATORY HOT PLATE

Our large expertise and market understanding help us manufacture and export superior quality range of Round Electric Hot Plate.

## About:

Manufactured in accordance with predefined industry norms with the help of latest technology by utilizing quality tested components, this hot plate is available in varied specifications as per the exact needs of the clients. Ideal for scientific researches, offered hot plate is used for providing heat to glassware and its content.



# Specification

## Of Laboratory Hot Plate



Laboratory hot plates are available in round and rectangular shapes.

Body is fabricated out of thick mild steel sheet duly finished in powder coated paint.

The hot plate made of cast iron, duly finished in heat resistant black paint is firmly mounted on the body.

Heavy duty heating elements are securely laid under the plate.

To operate on 220 V AC 50 Hz single phase.

Temperature is controlled by an energy regulator.

Provided with an indicator lamp for the mains and cord & plug.

## CONTACT:



Plot No. 6 - Ekta Vihar, Ambala Cantt.133001



Mobile No.: 8295732110, 7056732110



Sales@labsolutionindia.com



Labsolutionindia.com



**GeM**  
Government  
e Marketplace  
OEM Regd



**MSME**  
MICRO, SMALL & MEDIUM ENTERPRISES  
सूक्ष्म, लघु एवं मध्यम उद्यम