

Refrigerated Water Bath

Low Refrigerated Circulators for heating and cooling. The units are suitable for external temperature applications with simultaneous operation in the bath tank.

About:

- The units provide particularly high heating and cooling capacities for rapid heat-up and cool-down times even for large volume, external systems. With handle or castors to allow easy transport.
- With drain port for conveniently draining the bath fluid. The units offer an improved insulation to prevent ice-formation and have a filling level display.



Key Features

Of Refrigerated Water Bath



Construction	Inner SS and outer MS powder coated
Temperature range	0°-90°
Accuracy	5 degrees
Sensor	PT 100
Heaters	'L' type immersion heater
Cooling	Through hermetically sealed compressor
Refrigerant	R-134A
Water Circulation	Through pump
Power	220 volts AC

CONTACT:

-  Plot No. 6 - Ekta Vihar, Ambala Cantt.133001
-  Mobile No.: 8295732110, 7056732110
-  Sales@labsolutionindia.com
-  Labsolutionindia.com



GeM
Government
e Marketplace
OEM Regd

