

# ROTA MANTLE

The Rota mantle is a combination of heating Mantle cum Magn-etic Stirrer to obtain heating & stirring in flat bottom round flasks. The unit has two tier construction with the upper comp-artment consisting heating mantle made of fiberglass yarn & fitted with heating element.

## SPECIFICATION:

- Being a renowned firm, we are active-ly engaged in offering a wide range of Rota Mantles to fulfill the specific demands of our patrons.
- These offered products are exception-ally rigid in state and ensure longer operational life.
- All these products are precisely accur-ate and are compact in design.
- These products are highly resourceful and ensure smooth functioning.



# Key Features Of Rota Mantle



| MODELS      | Max. Speed<br>R.P.M. | Stirring<br>Capacity<br>Ltrs. | Heating<br>Capacity<br>Watts | Dimension<br>(W x D x H)<br>mm | Stirring<br>Paddle<br>f x L mm |
|-------------|----------------------|-------------------------------|------------------------------|--------------------------------|--------------------------------|
| LSI-RM -100 | 1200                 | 100 ml                        | 60                           | 140 x 140 x150                 | 7 x 15                         |
| LSI-RM -250 | 1200                 | 250 ml                        | 120                          | 150 x 150 x 160                | 7 x 15                         |
| LSI-RM- 500 | 1200                 | 500 ml                        | 200                          | 160 x 160 x 170                | 9 x 25                         |
| LSI-RM- 501 | 1200                 | 1 ltr                         | 200                          | 190 x 190 x 190                | 9 x 25                         |
| LSI-RM- 502 | 1200                 | 2 ltrs.                       | 400                          | 230 x 230 x 230                | 9 x 35                         |
| LSI-RM- 503 | 1200                 | 3 ltrs.                       | 500                          | 250 x 250 x 250                | 9 x 35                         |
| LSI-RM- 505 | 1200                 | 5 ltrs.                       | 2 x 300                      | 290 x 290 x 290                | 13 x 50                        |
| LSI-RM- 510 | 1200                 | 10 ltrs.                      | 2 x 500                      | 400 x 400 x 340                | 13 x 75                        |
| LSI-RM- 520 | 1200                 | 20 ltrs.                      | 3 x 500                      | 450 x 450 x 380                | 13 x 75                        |

## CONTACT:

 Plot No. 6 - Ekta Vihar, Ambala Cantt.133001

 Mobile No.: 8295732110, 7056732110

 Sales@labsolutionindia.com

 Labsolutionindia.com



**GeM**  
Government  
e Marketplace  
OEM Regd



**MSME**  
MICRO, SMALL & MEDIUM ENTERPRISES  
सूक्ष्म, लघु एवं मध्यम उद्यम