

Wrist Action Shaker

It has a heavy cast iron base with rugged construction and two side arms for holding flasks. Fitted with the heavy-duty motor with variable speed can hold flask from 100 to 1000ml. Capacity to Hold 2 Flask, 4 Flask or 8 flask.

About:

We are engaged in offering wide range of Wrist Action Shaker. Our experienced designers have used their vast knowledge and advance technology to make it a worthy product at market leading price. Client can avail these from at industry leading.

- Resembles hand mixing and is an appropriate solution for dissolution, extractions & homogenization.
- Easily Detachable arms.
- Timer option for accurate mixing time.
- Digital RPM controller available for exact results.



Key Features

Of Wrist Action Shaker



Body	Made of Mild Steel Powder coated
Motor	Heavy duty geared DC motor
Speed	Variable Speed
Shaking Angle	0 to 15°
Speed	20 to 450 RPM
Max Load	10 lb.
Clamps	2 Clamps to hold four Vessels
Oscillations	Max, 450 per minute
Timer	1 to 60 Minutes Timer or Continuous running
Unit Dimensions	17.5" L x 17.5" W x 12" H
Unit Weight	@ 15 Kgs.
Mounted on	Heavy duty Vacuum Cups
Power Supply	220 V Single Phase

CONTACT:

-  Plot No. 6 - Ekta Vihar, Ambala Cantt.133001
-  Mobile No.: 8295732110, 7056732110
-  Sales@labsolutionindia.com
-  Labsolutionindia.com



GeM
Government
e Marketplace
OEM Regd

