

# ANALOG TYPE OIL BATH

These oil baths can be availed with added accessories such as GMP model, and digital indicator and controller for constant temperature. Further, it is provided with a heavy duty motor that makes its operation swift and easy.

## About:

### Features:

- Long working life.
- Excellent strength.
- Corrosion resistance.



# Key Features

## Of Analog Type Oil Bath



Double walled: inner S.S. And outer made of M.S. Powder coated and consists of welded inner tank made of stainless steel 304.

Immersion type heater.

Temperature controlled by thermostat.

Options: GMP model, and digital indicator and controller for constant temperature.

Rectangular model with oil circulating motor.

These oil baths can be availed with added accessories such as GMP model, and digital indicator and controller for constant temperature.

## CONTACT:

 Plot No. 6 - Ekta Vihar, Ambala Cantt.133001

 Mobile No.: 8295732110, 7056732110

 Sales@labsolutionindia.com

 Labsolutionindia.com



**GeM**  
Government  
e Marketplace  
OEM Regd

