



CAMERA LUCIDA MIRROR TYPE

CAMERA LUCIDA MIRROR TYPE

- ✔ Easily trace the image of the microscope object on the plain paper
- ✔ The artist sees both scene and drawing surface simultaneously
- ✔ Aiding in the accurate rendering of perspective
- ✔ Can be used by school, college students etc
- ✔ Package content: 1 mirror type lucida camera

About Camera Lucida Mirror Type

Mirror type, to trace the image of a microscopic object on a paper. A beam splitter prism of optical glass is positioned over the eyepiece with a large size swinging mirror which causes the image of the object to appear on a plane surface enabling its tracing.

