



# **DISSECTING MICROSCOPE**

# DISSECTING MICROSCOPE

- ✓ Base Heavy round base with precisely designed body.
- ✓ Focussing Sensitive focussing, done by rack and pinion arrangement.
- ✓ Revolving arm provided for moving magnifying lens overfull stage area.
- ✓ Stage 85mm x 75mm with glass plate.
- ✓ Illumination A plano-concave lens fitted in fork for light reflection.
- ✓ Optics 10x & 20x Eyepieces

## About Dissecting Microscope

Dissection Microscope, Heavy, round C.I. base. All brass C.P. parts like pipe, knobs, etc. Eye-pieces 10x & 20x are highly polished and made out of optical glass & an easily replaceable 83mm Bull-lens (most -ly used for examination of insects). Rack & pinion motion for focusing. Mirror illumination. In an attractive & thick plywood box.

