



# **JUNIOR MICROSCOPE** **MOVING CONDENSOR**

# JUNIOR MICROSCOPE MOVING CONDENSOR

- ✓ Stable & Robust, traditional Horse shoe base.
- ✓ Huygenian eyepiece 10x or 15x.
- ✓ Body inclinable through 90°, Tube Length 160mm.
- ✓ Triple nosepiece.
- ✓ A1 series chromatic objectives 10x, 40x (SL).
- ✓ Fixed Stage size 110 x 110mm with stage clips to hold the slide.
- ✓ Moving condenser by Rack and pinion, condenser N.A 1.25, aperture selection through a Diaphragm.
- ✓ Seperate coarse and fine focusing knobs. Coarse movement by Rack & Pinion and fine movement by screw lever mechanism.
- ✓ Illumination system Plano concave reflector mirror mounted on gimbal mount.

## About Junior Microscope Moving Condensor

### Optional Accessories

Wooden Box

2/3 Megapixel Digital Camera for Microscope

Pin Type Lamp (Bulb/Led) 220v

