

LABORATORY REFRIGERATOR

We are affianced in offering Laboratory Refrigerators with counter height and high capacity. Offered refrigerators are broadly acknowledged by our patrons as they provide utmost storage space for reagents, vaccines and other critical samples. Sturdy construction, proper cooling and long life are some of the unique attributes.

These laboratory refrigerators are chosen worldwide for their reliability, accurate temperature control, security benefits and usability. As expected from Bionics Scientific a leading laboratory refrigerator manufacturer in India, the refrigeration system is environment friendly and only CFC free compressors that are virtually noiseless in nature used in these systems. High grade PUF insulation minimizes heat loss.

Key Features :

- Inner body made up of medical grade compression molded plastic
- Automatic defrost for frost free operation
- Uniform temperature with forced air circulation
- Clear product visibility with glass door
- Adjustable coated wire mesh shelves
- A range of sizes are available from 50 litres to 1000 litres or more with either single or double outer doors.



Specification

Of Laboratory Refrigerator



Capacity	400 Liters
Temperature Range	2 to 8 Degree Celsius
Product Dimensions (WDH)	610 x 584 x 1981
Number of Trays	4
Door Type	Vertical Door
Controller	PID Controller
Refrigerant	HFC (CFC-Free)
Construction	Exterior Powder Coated Mild Steel / SS (Optional) Interior Stainless Steel
Insulation	PUF insulation
Power Supply	230 V, 50 Hz

CONTACT:

-  Plot No. 6 - Ekta Vihar, Ambala Cantt.133001
-  Mobile No.: 8295732110, 7056732110
-  Sales@labsolutionindia.com
-  Labsolutionindia.com



GeM
Government
e Marketplace
OEM Regd

