



MONOCULAR MICROSCOPE

MONOCULAR MICROSCOPE

- ✔ Frame– All aluminum pressure die-casted body.
- ✔ Paint– Powder coated, available in various desired shades.
- ✔ Nosepiece– Quadruple ultra smooth nosepiece.
- ✔ Focusing system– Separate knobs for coars and fine motion having the least count of .002 mm.
- ✔ Stage– Cxl mechanical stage works on dovetail slides.(other stages also available).
- ✔ Condenser Bright field sub-stage condenser having 1.2 N.A.
- ✔ Eyepieces & objectives– Set of 4x, 10x, 40x, and 100x and a paired set of 10x wide field eyepieces.
- ✔ Illumination– 6v,20w halogen lamp needs 220v and 50Hz supply.
- ✔ Instruction Manual, Dust cloth, Dust cover.
- ✔ standard configurations
- ✔ Eyepieces & objectives- Set of 4x, 10x, 40x, and 100x and a paired set of 10x wide field eyepieces.

About Monocular Microscope

Optional accessories

- Wooden Box
- 2/3 Mega Pixel Digital Camera for Microscope

