

# PLATELET AGITATOR

Platelet Agitator model is specifically designed to maintain donor platelets an even suspension throughout the blood plasma. Platelet agitator agitates at the fixed speed of 70 to 80 strokes per minute.

01 > 02 > 03 > 04

Platelet Agitator agitates at the fixed speed of 70 to 80 strokes per minute.

The self are pure in stainless steel.

Maximum 48 standard blood bags can be accommodated.

Fixed speed agitation with optimum amplitude for better viability of platelets.



# Specification Of Platelet Agitator



Brand	Lab solution India
Fixed Speed	70 to 80 strokes per minute
Number of Trays	9
Material	Stainless Steel
Size	38 x 45 x 40 cm
Voltage	220 V
Freequency	50 Hz

## CONTACT:

 Plot No. 6 - Ekta Vihar, Ambala Cantt.133001

 Mobile No.: 8295732110, 7056732110

 Sales@labsolutionindia.com

 Labsolutionindia.com



**GeM**  
Government  
e Marketplace  
OEM Regd

