



# PROJECTION MICROSCOPE

# PROJECTION MICROSCOPE

- ✓ Monocular Fixed Straight Tube, length 160mm.
- ✓ Nosepiece Triple revolving nosepiece with accurate centering and positive click stops.
- ✓ Stage Mechanical stage of 120x120mm.
- ✓ Illumination 20W lamp, 220/110Volts. Lamp slots provided for natural cooling of the illuminator.
- ✓ Objectives Achromatic 10X-40XSL.
- ✓ Eyepieces Huygenian 10X
- ✓ Screen 150mm diameter Glass Projection Screen is provided, can be rotated over 360-degrees on Microscope tube.
- ✓ Included Items
  - Includes: Instruction Manual, Dust cloth, Dust cover,
- ✓ standard configurations
  - Objectives Achromatic 10X-40XSL.
  - Eyepieces Huygenian 10X

## About Projection Microscope

### Optional accessories

- Wooden Box
- Mechanical Stage

