

RECIPROCATING SHAKING MACHINE

Reciprocating shakers by Lab Solution India are ideal for mixing blood samples, chemical extractions and silver staining of proteins in polyacrylamide gels. Also called Kahn Shaker or reciprocating mixer machine, these units are fitted with analog control system; variable speed control enables you to operate the machine with gentle or vigorous mixing speed. In addition; these reciprocating shakers can also be used in warm rooms, environmental chambers, large refrigeration units and incubation systems. Platform is adjustable, provides you flexibility to change dimensions of compartments by adjusting rods. Our reciprocating shakers come with number of optional accessories in order to make the machine suitable with variety of applications; different sizes of flasks, bottles, tube racks are available on request.

About:

Durability is one of the most exciting features that stand out our reciprocating shakers among other models available on the Indian market today. These units are built for years of long services life; power coated Mild Steel (MS) construction ensures durable performance anywhere they used. To know more about our reciprocating shaker models including price, scope of customization and availability, please submit your query, our sales executives will assist you with all the information you need.



Key Features

Of Reciprocating Shaking Machine



Speed	50 to 250 rpm
Motion	Reciprocating
Control	Variable speed control
Construction	Powder Coated MS
Power	220 / 230 Volts 50 Hz
Optional	- 0 - 60 minutes mechanical timer - Digital timer - Digital rpm - Tube racks - Bottle holder

Reciprocating flow motion	Adjustable Kahn racks
Variable speed control	Long service life
Available in choices of platforms	Heavy duty steel construction

CONTACT:

 Plot No. 6 - Ekta Vihar, Ambala Cantt.133001

 Mobile No.: 8295732110, 7056732110

 Sales@labsolutionindia.com

 Labsolutionindia.com



GeM
Government
e Marketplace
OEM Regd

