

RECTANGULAR WATER BATH

12 Holes (Digital Model)

The electrically operated Rectangular Water Bath are made of thick stainless steel and feature constant water levels maintained inside them. For optimum performance, these also come complete with concentric rings having 75 mm dia, wire plugs and with suitable drainage facility. The bath machine also comprises thermostatic controller that allows better control of temperature levels with removable glass tray and stainless steel platform for drying of sections. Some of the applications these are used in include for sterilizing purposes. These water baths also come fitted with heating element for uniform heating throughout the system and finds usage in areas like medical institutions, research laboratories, agricultural institutions, bio-chemical industries as well as in many other industries.

Product Description

The general purpose digital Rectangular Water Bath has a temperature regulated sensor and is essential for

- Drying,
- Condensing
- Distilling
- Soaking chemical Agents
- Drug
- Biological Products Extensively.

It can also be used for heating and other temperature related experiments.



Specification

Of Rectangular Water Bath 12 Holes
(Digital Model)



Voltage	220 / 230 V
Temperature Range	Ambient to 100 Deg C
Capacity	12 holes
Heating Power	1000 W
Temperature Accuracy	+/- 1 Deg C
Inner Wall	Stainless steel 304 grade
Outer Wall	Mild Steel
Top	Stainless Steel
Insulated By	Glass wool to avoid thermal loss
Depth of Wells	4.7 inches
Temperature Control	By Thermostat
Minimum Order Quantity	1

CONTACT:

 Plot No. 6 - Ekta Vihar, Ambala Cantt.133001

 Mobile No.: 8295732110, 7056732110

 Sales@labsolutionindia.com

 Labsolutionindia.com



GeM
Government
e Marketplace
OEM Regd

