



# **SINGLE NOSE MICROSCOPE**

# SINGLE NOSE MICROSCOPE

- ✓ Stable & Robust, traditional Horse shoe base
- ✓ Huygenian eyepiece 10x or 15x (any ONE )
- ✓ Body inclinable through 90°, Tube Length 160mm
- ✓ A1 series Achromatic objectives 10x
- ✓ Fixed Stage with stage clips to hold the slide
- ✓ SMOOTH coarse focusing knobs. Coarse movement by Rack & Pinion
- ✓ Illumination system Plano concave reflector mirror mounted on gimbal mount

## About Single Nose Microscope

### Included Items

Includes: Instruction Manual, Dust cloth, Dust cover.  
**standard configurations**

- A1 series chromatic objectives 10x.
- Huygenian eyepiece 10x or 15x (any ONE ) Optional accessories.
- Eye Piece.
- Light Source.

