

VDRL SHAKER

(Digital Model)

VDRL shaker is an essential equipment for a Venereal Disease Research Laboratory, where it is used for the mixing of syphilis testing card, emulsion reaction, and complement fixation reaction while can also be used to blends of clinical diagnostic tests requiring gentle agitation. Bionics Scientific, designs and manufactures VDRL Shakers that are compact, durable and offer fixed as well as variable speed controls and can accommodate wide range of VDRL Slides, Small Flasks, Bottles & Beakers.

Features:

- Compact design, durable construction
- Smooth mixing
- Platform with holder clamp and springs
- Heavy duty motor
- Highly balanced drive mechanism
- Non-slip rotating base
- Ideal for serological & syphilis tests
Uniform speed



Specifications

Of VDRL Shaker (Digital)



Shaking speed 40 to 200 rpm

Amplitude 10 mm

Timer can be set up to 999 minutes

Micro controller based programme

Brushless Motor

Interchangeable Platform

LED display

CONTACT:

 Plot No. 6 - Ekta Vihar, Ambala Cantt.133001

 Mobile No.: 8295732110, 7056732110

 Sales@labsolutionindia.com

 Labsolutionindia.com



GeM
Government
e Marketplace
OEM Regd

