

BLOOD DONOR CHAIR

Solid surface seat to ensure maximum comfort for the donor. The seat has a premium rated foam filled back rest supported by a steel structure. The sturdy construction of this seating ensures that the chair can be used year after year

About:

Comfortable chair type with soft padding for cushioning and rexin cover. Seat, back rest and leg rest size designed for donor comfort. The Blood Donor chairs is used in various application such as Blood collection chair for collection of Blood, Blood donation chair for donation of Blood.



Key Features

Of Blood Donor Chair



Framework made of Stainless Steel.

Back rest & leg rest section adjustable on Ratchet / Lever.

Legs fitted with PVC Stumps.

5 cms foam cushioned top covered with rexine.

Head section & leg section move simultaneously.

Adjustable Arm Rest.

CONTACT:

 Plot No. 6 - Ekta Vihar, Ambala Cantt.133001

 Mobile No.: 8295732110, 7056732110

 Sales@labsolutionindia.com

 Labsolutionindia.com

