

CABINET DRYER

Cabinet Dryer consist of an enclosed chamber where products to be dried are placed on stainless steel trays. These drying systems reduce the relativity humidity allowing the moisture to evaporate quickly. Cabinet dryers are widely used for redu-
-cing the moisture content of sensitive products.

01 > 02 > 03 > 04

Laboratory Dryers are constructed of powder coated mild steel.

Each shelves are removable and adjustable.

Drying temperature is set by a thermostat. The thermal safety protection with an independent high temperature.

Limit cut-off thermostat can protect the samples from destruction.



Key Features Of Cabinet Dryer



Capacity	400 ltr
Shelves	5 shelves
Temperature	30 to 70 degree
Tray Material	Stainless Steel
Outer Material	Mild Steel
Air-Flow Direction	Horizontal
Humidity	10 percent
voltage	220/380/460 , 50/60Hz

CONTACT:

-  Plot No. 6 - Ekta Vihar, Ambala Cantt.133001
-  Mobile No.: 8295732110, 7056732110
-  Sales@labsolutionindia.com
-  Labsolutionindia.com



GeM
Government
e Marketplace
OEM Regd

