

Portable Autoclave (Aluminium)

- All body made of Aluminium
- Ring and Lid made of aluminium Casting.
- Used for sterilization Surgical instruments and dressing.
- Material – Aluminium Hydraulically tested upto 30 psi.
- Fitted with pressure gauge.
- The unit has spring loaded safety valve of pressure adjustment and a dead weight safety valve

INCLUDING:

Material	Aluminium
Orientation	Vertical
Pack Type	Box
Usage	Laboratory
Capacity	20 L, 25L



Technical Specification of Portable Autoclave (Aluminium)



Model	Size(mm)	Size(inches)	Capacity(liter)
901	300X300	12X12	20ltr
901	300X375	12X15	25ltr

CONTACT:



Plot No. 6 - Ekta Vihar, Ambala Cantt.133001



Mobile No.: 8295732110, 7056732110



Sales@labsolutionindia.com



Labsolutionindia.com



GeM
Government
e Marketplace
OEM Regd



MSME
MICRO, SMALL & MEDIUM ENTERPRISES
सूक्ष्म, लघु एवं मध्यम उद्यम