

UNIVERSAL WATER BATH

Universal Water Bath feature double walled construction with outer chambers designed using duly powder coated mild steel finish as well as inner bath chamber and lids design -ed using superior grade stainless steel so as to provide for optimum functional support in involved procedures. The filling of spaces between walls using mineral wool support also helps in avoiding thermal losses. Further, the presence of temperature controls delivering accuracy of $\pm 0.5^{\circ}\text{C}$ as well as system's availability with 6 holes/ 12 holes of 75mm concentric rings also makes these suitable for diverse lab procedures.

Key Features :

- Displays actual and set point temperature simultaneously.
- Integrated timer.
- Reservoir sizes from 2 to 28 liters.
- Reservoir made of 304 Stainless Steel.
- Reservoir drain (10 liter and larger baths).
- Primary and automatic safety thermostats.
- Also compatible with waterless bath beads.



Specification

Of Universal Water Bath



Capacity	10 Litre
Material	SS
Voltage	230 V
Temperature Range °C	Ambient +5° to 100°
Phase	Single Phase

STANDARD SIZES AVAILABLE

Size Inside (SS)	300 x 225 x 175mm	450 x 300 x 175mm
Size In Inches	12 x 9 x 7	181x 12 x 7
Capacity in Liters	14 Liters	25 Liters

CONTACT:

 Plot No. 6 - Ekta Vihar, Ambala Cantt.133001

 Mobile No.: 8295732110, 7056732110

 Sales@labsolutionindia.com

 Labsolutionindia.com



GeM
Government
e Marketplace
OEM Regd

