

AUTOCLAVE/STERILIZERS

AUTOCLAVE TABLE TOP

CONSTRUCTION AND SALIENT FEATURES

- Bench-top, pre and post vacuum autoclave with automatic temperature control system, self-test system before operation, self-test mode diagnose, fully automatic steam ejection and safety device containing safety valve, automatic sterilization cycles and auto detections for steam sterilization of glass wares and plastic wares.
- Tank Capacity: 20 to 25 liters.
- Stainless steel seamless chamber and epoxy coated chamber cover for more safety.
- Sterilize cycle, dry timer and indicator lamps display on panel which allows the operator to visibly acknowledge the phase. It is having programmable based controller.
- Safety locking mechanism that prevents opening while the autoclave is under Pressure. Complete control panel is provided in the front.
- Steam release valve will exclude cool air and extra steam.
- If electrical AC power is unavailable or is interrupted, the autoclave shall automatically shut off power and buzz to protect autoclave and instruments
- The autoclave controls are easy to operate. Just choose the mold and time you want first and then press start.
- It is thoroughly tested and proven technology. (CE 1023 approved).
- Uniform rubber gasket which endure high pressure and temperature with low dilatibility.



SPECIFICATIONS

- Sterilize temperature: 121°C or 134°C
- Electrical: AC 220V or 110V, 2100 watts / 220V and 1600 watt / 110 V
- Water reservoir capacity: Approx. 5 L
- Safety device: safety valve, pressure regulator, and overheat protection
- Maximum Temperature: 150°C
- Maximum Pressure: 0.25mpa
- Warranty: One year