

FUME HOOD

Used when performing acid digestion applications. The fume hood is designed in order to throw out toxic or harmful Fumes. Fume hood are used to protect the persons and eliminate the toxic gases through ducting by exhaust blower.

01 > 02 > 03 > 04

Easy-of-Use: Sliding sash movement offers easy and thorough cleaning to every users.

The high-grade transparent glass door offers excellent visibility of the chemical handling.

The most suitable UV light provides you with optimum decontamination in the work zone surface.

Highly effective backup safety device for the contaminant and exhaust of toxic, offensive or flammable materials.



Features Of Fume Hood



Power	2hp, 3 Phase
Dimension	15mm
Pressure	2.5 W.G.
Inner Size	2'x2'x2' or 3'x2'x2' or 4'x2'x2'
Front Door Material	6mm Acrylic Sheet
Volume	1800 CFM
Motor Speed	1440 RPM
Type	Centrifugal
MOC	1 mm Mild Steel

CONTACT:

-  Plot No. 6 - Ekta Vihar, Ambala Cantt.133001
-  Mobile No.: 8295732110, 7056732110
-  Sales@labsolutionindia.com
-  Labsolutionindia.com



GeM
Government
e Marketplace
OEM Regd

